

One Man Sling

Manual English

 etac®



Illustrations

Illustration 1





Illustration 2



Contents

	Heading	Page
1.	General/Features.....	4
2.	Practical handling.....	5
3.	Handling/Materials/Cleaning	6
4.	Size/Combinations	7

General

Thank you for choosing an Etac product.

To avoid accidents and injury when moving and handling the products, read the manual carefully.

The person described as "the client" in this user manual is the person who is lying or sitting on the product. The carer or carers are the people who manoeuvre the product.



This symbol appears alongside the text in the manual. It draws the reader's attention to points at which there may be a risk to the health and safety of the client or carer.

The products comply with the standards applicable for Class 1 products in the European Council Directive MDD 93/42/EEC on medical devices.

At Etac we strive to improve our products all the time and therefore we reserve the right to make changes to products without prior warning. All measurements given on illustrations and similar material are for guidance only and Etac cannot be held liable for errors and defects.

The information given in this manual, including recommendations, combinations and sizing, does not apply to special orders and modifications. If the customer makes adjustments, repairs or combinations not predetermined by Etac, the Etac CE certification and Etac warranty will be voided. If in doubt, please contact Etac.

Warranty: Two-year warranty on material and manufacturing defects, provided that the product is used correctly.

For further information on the Etac transfer range, see www.etac.com.

The product can be scrapped in accordance with national regulations.



Washing



Iron



Never leave on the floor



Max. weight



Handwash



Do not dryclean



Warning



The product can be scrapped in accordance with national regulations



Wipe off



Do not bleach



Read the User Manual



Manufacturer



Tumbledry



Risk of sliding down



Batch no./Lot no.



Class 1 - European Council Directive MDD 93/42/EEC on medical devices

Features

The One Man Sling has been specially developed for manual handling a client higher up in the bed when only one carer is available for the task.

Use

Positioning

The black webbing belt is fastened to a stable part of the bed frame at the head end.

Attach the sheath to the webbing belt and feed it past the handles, so that it glides without obstruction.

Use

Moving the client higher up in the bed

The client bends his/her legs slightly.

The One Man Sling is pushed under the client's legs and is pulled until the sheath lies against the client's rear.

The rest of the webbing strap with the handles is fed to the opposite side of the bed.

The carer moves around to the opposite side of the bed at the head end. The carer positions a Multiglide under the client's pillow and pulls down so that it covers the client's shoulders.

Further glide devices can be added under the client's rear (and under the legs if the client is not able to bend at the knee).

The carer stands with one foot forward, grasps one or two of the handles closest to the client's hip and performs weight transference (Illustrations 1 and 2). The procedure can be repeated if necessary.

After manual handling, remove the Multiglide and feed the One Man Sling back to the opposite side at the head end, where it can be stored until the next time.

Legs into bed

Attach the One Man Sling to a stable part of the bed frame at the foot end on one of the long sides (the same side as where the client sits on the edge of the bed).

Feed the One Man Sling under the client's knees (closest to the foot end) and around the ankle of the other leg, so that the client's legs are above the One Man Sling.

When the client is ready, the carer pulls the handles of the One Man Sling in the direction of the foot end.

Both the client's legs are now pulled towards the foot end. The carer continues to move around the end of the bed in an arch, ending up at the opposite side of the bed.

There is no need to lift the webbing belt. The carer simply has to move in an arch around the foot end.

When the carer reaches the foot end on the opposite side, she/he can continue to pull on the strap until the client's legs are in the middle of the bed.

For optimal gliding movement, a Glide Cushion or Multiglide can be used at the edge of the bed at the foot end under the One Man Sling.

NB! If the client is unable to lay down the upper body by him/herself, this manual handling must be carried out by two carers. The bedhead can be raised if necessary.

Use

Special features

Always check the product before use and after washing.

Never use a defective product. If the product shows signs of wear-and-tear, it must be scrapped.

Read these instructions carefully.

It is important that carers receive instruction in manual handling.

Etac offers advice and training for carers.

For further information, contact Etac.

Make sure you know exactly how much weight you are trying to move.

Call for extra staff if necessary!

Always use the correct manual handling techniques.

Remember to use the client's own resources.

Never leave One Man Sling on the floor.



To ensure that the client feels safe and that every manual handling is smooth, always plan the manual handling in advance.

If in any doubt - please contact Etac.

Materials

Belt: 100% polyester

Sheath:

Inner surface: 100% nylon.

Fill: 100% polyester (machine-washable at 95°C)

Outer surface: polyurethane

Cleaning

The product should be checked regularly, preferably each time it is used.

Check that seams and fabric are not damaged.




Never try to repair a product yourself. Never use a defective product.

Washing instructions:

Separate the two parts before washing.

Do not use fabric softener.



Articles

Article no.	Description	Size (cm)
IM460	One Man Sling, 6 handles	W7.5 x L300

Combinations

Article no.	Description	Size (cm)
IM460	One Man Sling, 6 handles	W7.5 x L300
IM40	Multiglide, waterproof, blue strap	B40 x L40
IM66	Multiglide, waterproof, yellow strap	B60 x L60
IM77	Multiglide, waterproof, brown strap	B72 x L70
IM80	Multiglide, waterproof, red strap	B72 x L100
IM82	Multiglide, waterproof, light brown strap	B70 x L200
IM83	Multiglide, waterproof, green strap	B72 x L130
IM88	Multiglide, waterproof, green strap	B80 x L80
IM0077	Multiglide, black strap	B72 x L70
IM0080	Multiglide, black strap	B72 x L100
IM0083	Multiglide, black strap	B72 x L130
IM0088	Multiglide, black strap	B80 x L80
IM60/50	GlideCushion, green, cotton/polyester	B60 x L50
IM60/50N	GlideCushion, green, nylon top side	B60 x L50
IM60/90	GlideCushion, green, cotton/polyester	B60 x L90
IM69/60	GlideCushion, green, cotton/polyester	B69 x L60
IM69/70	GlideCushion, green, cotton/polyester	B69 x L70

Etac Sverige AB
Box 203, 334 24 Anderstorp, Sweden
Tel +46 371-58 73 00
Fax +46 371-58 73 90
info@etac.se www.etac.se

Etac AB (export)
Box 203, 334 24 Anderstorp, Sweden
Tel +46 371-58 73 30
Fax +46 371-58 73 90
info@etac.se www.etac.com/exp

Etac AS
Pb 249, 1501 Moss, Norway
Tel +47 815 69 469
Fax +47 69 27 09 11
hovedkontor@etac.no www.etac.no

Etac A/S
Egeskovvej 12, 8700 Horsens, Denmark
Tel +45 79 68 58 33
Fax +45 75 68 58 40
info@etac.dk www.etac.dk

Etac GMBH
Bahnhofstraße 131, 45770 Marl, Germany
Tel +49 236 598 710
Fax +49 236 598 6115
info@etac.de www.etac.de

Etac Holland BV
Fluorietweg 16a, 1812RR Alkmaar, The Netherlands
Tel +31 72 547 04 39
Fax +31 72 547 13 05
info.holland@etac.com www.etac.com